

**UNDERSTANDING LOW VISION TELESCOPES**

*Basic and Advanced*

**Ian L Bailey** *OD, MS, DSc, FCOptom, FAAO*  
 School of Optometry  
 University of California, Berkeley

*ibailey@berkeley.edu*

**Applications for Telescopes in Low Vision**

**Observation tasks from audience**

- Theater
- Ballet
- Movies
- Sports
- Classroom, Public speakers
- TV in communal environments

**General tasks**

- Scenery
- Faces
- Cars
- McDonalds Menu
- Computer screens
- Top shelf items

**Mobility tasks**

- Traffic signals
- Street signs
- Bus numbers, airport signs
- Building directories, wall signs
- Checking traffic

**Choices For the Clinician**  
*or the Patient*

**Optical Factors**

- Magnification
- Field size
- Exit pupil size or entrance pupil for stars
- Adjustable focus/fixed focus
- Range of focusing (for near and ametropia)
  - Can you include astigmatic Rx ?
  - Image quality
  - Lens Coatings

**Design Features**

- Binocular/Monocular
- Hand-held /Head mounted
- If head mounted Bi-optic ?
- Tripod/monopod
- Eye-cup, Objective shroud
- Neck, wrist strap, handle
- Influenced by frequency of use*  
*(Protracted, Intermittent, Occasional)*

**Practical Considerations**

- Cost
- Durability
- Portability (*weight/bulk/holder*)
- Ease of use (*mechanical ergonomics*)
- Cosmesis (*conspicuity and perceived acceptability*)

**Course Outline**

**Review of basic theory**

- Afocal telescopes
- Focal telescopes *for ametropia or close focus*

**Apertures and Stops**

- Image brightness
- Field of View *and Field of Fixation*

**Vergence amplification**

- Practical implications
- Checking magnification
- Measuring BVP

**Fitting bioptic telescopes**

**Emerging technologies** *zoom, autofocus, video*

**How do afocal telescopes work?**  
 Parallel rays in gives parallel rays out

**Objective**


A converging element (Convex lens, Concave mirror) forms a real inverted image.



*Weak objective- long focal length - larger inverted image*

**Ocular or Eyepiece**

Objective 's image must be at primary focal plane of ocular

*Stronger ocular - short focal length - large image (infinity)*





### FOUR BASIC TELESCOPE FORMULAE

The two obvious relationships

$M_{ts} = - P_{oc} / P_{obj}$  Ratio Ocular / Objective

$d = (1 / P_{obj}) + (1 / P_{oc})$  Sum of focal lengths


Two similar formulae (similar to spectacle magnification formula)

$M_{ts} = 1 / (1 - d P_{obj})$

$M_{ts} = 1 - d P_{oc}$

### Comparing Galilean and Keplerian

	GALILEAN	KEPLERIAN
Optical path length	Short	Long
Image	Erect (M=+ve)	Inverted (M= -ve)
Prisms	No	Yes <i>Makes it terrestrial</i>
Ocular	Single lens	2 or more lenses
Optical surfaces	4	10 or more
Exit pupil	<i>inside</i>	<i>outside</i>
Field of view	<i>smaller</i>	<i>larger</i>
Field limits	<i>tapered</i>	<i>sharp outline</i>
Image quality	<i>poorer</i>	<i>better</i>



### AFOCAL TELESCOPES

AFOCAL if Parallel rays in  
give Parallel Rays out

*Afocal telescopes*  
have some mathematically sweet properties

*Usually telescopes not quite afocal*  
Objects not at infinity  
Patients often do not want or demand image at infinity

### Adapting Afocal Telescopes for Ametropia

*Assume object at infinity*  
*What can you change?*

- (a) **OCULAR** Place ocular against spec Rx  
OR Change ocular to include spec Rx
- (b) **LENGTH** Change length, changes exit vergence
- (c) **OBJECTIVE** Low Powered Lens Cap OR power change at objective changes exit vergence **RARELY USED**

### (a) POWER CHANGE AT OCULAR

*Simple and intuitive*

Hold afocal telescope against own spectacle lens  
OR  
Change ocular lens by adding spectacle Rx

Magnification gain = magnification of original telescope

**Total magnification = (Spectacle mag) (Telescope mag)**

No differences for *Galilean/Keplerian*

**EVD = viewing distance /  $M_{ts}$**

### (b) CHANGE LENGTH

Changing length, changes exiting vergence


Increasing length adds plus (**HYPEROPE**)  
*Also increase length to focus for near (adding plus)*

Decreasing length adds minus (**MYOPE**)

### FOR HYPEROPE Increase length

*Re-assign PLUS power from ocular to provide Rx.*

**Keplerian telescope**  
Borrow plus from ocular  
*Weaker ocular, longer TS*  
**Reduced magnification**




**Galilean telescope**  
Borrow plus from ocular  
*Stronger ocular, longer TS*  
**Increased magnification**

*Be able to calculate TS length and  $M_{ts}$*

### FOR MYOPE Decrease length

*Re-assign MINUS power from ocular to provide Rx.*

**Keplerian telescope**  
Borrow minus from ocular  
*Stronger ocular, shorter TS*  
**Increased magnification**



**Galilean telescope**  
Borrow minus from ocular  
*Weaker ocular, shorter TS*  
**Reduced magnification**

*Be able to calculate TS length and  $M_{ts}$*

### TELESCOPES FOR NEAR VISION

- 1 Afocal TS with a cap  
$$EVD = u/M_{ts} = f_{cap}/M_{ts}$$
- 2 Adjust TS length of afocal telescope  
$$EVD = (u/M_{ts}) - f_{oc} \quad EVD \approx u/M_{ts}$$
  
Where  $f_{oc}$  = focal length of ocular (*usually 1-2 cm*)
- 3 Afocal telescope over a high add  
Rarely used  
$$EVD = M_{ts} \cdot f_{add}$$
  
Working distance =  $M^2 f_{add} - dM_b$   
System gives low magnification  
Long EVD and much longer WD

### Increasing length a 4x12 TS focused for 20 cm

Increase length to "add plus" for near vision focus.

4x12 Monocular telescope  $M = -100/25 = -4x$   
Tube length = 5 cm when afocal

Effectively borrow +5D from objective

Adjusted to focus on 20 cm  
Now system is like a 5x TS with a +5D cap  $M = -100/20 = -5x$

Adjusted to focus on 20 cm  
Now system is like a 5x TS with a +5D cap  $M = -100/20 = -5x$   
Light diverging from 20 cm  
Tube length now 6 cm  
EVD = 20 cm/5x = 4 cm    EVP = 25D

System behaves like an afocal telescope with a cap on the objective

### Labeling Telescopes

Two numbers  
Magnification x Objective diameter mm  
Manufacturers never state negative magnification for Keplerian TS's

4x12  
is a 4x TS, with objective 12 mm in diameter

10x50  
is a 10x TS with objective 50 mm in diameter

### EXIT PUPILS

Small images of Objective lens  
Located quite close to eyepiece lens  
Keplerian telescopes (Exit pupil outside TS)  
Galilean telescopes (Exit pupil outside TS)

Relevant to Image Brightness

Relevant to Field of View

Can help understand imaging by telescopes  
An important reference point

Can be used to determine magnification of TS

### EXIT PUPIL OF KEPLERIAN TELESCOPE

OBJECTIVE  $A_{obj} = \text{diam}$

EYEPIECE  $A_{eye} = \text{diam}$

EXIT PUPIL = image of objective

Exit pupil is a real image suspended in space  
Its size is  $A_{obj} / M_{TS}$   
It is located  $d / M_{TS}$  from the eyepiece lens

### EXIT PUPIL OF GALILEAN TELESCOPE

OBJECTIVE  $A_{obj} = \text{diam}$

EYEPIECE  $A_{eye} = \text{diam}$

EXIT PUPIL = image of objective

Exit pupil is a virtual image inside telescope  
Its size is  $A_{obj} / M_{TS}$   
It is located  $d / M_{TS}$  from the eyepiece lens

### Image Brightness and Telescopes

REFLECTIONS AND ABSORPTION

All telescopes reduce brightness of real objects  
Light loss (about 15-40%) from reflections from internal optical surfaces

Keplerian more light loss  
More surfaces (compound eyepiece and prisms)


Surface coatings reduce light loss

**BUT**  
TELESCOPES ENHANCE BRIGHTNESS OF POINT SOURCES (STARS)

### STARS OR POINT SOURCES

*Retinal image is smaller than a receptor diameter*


**For Point Sources**  
 As the eye gets closer to a point source, the pupil becomes larger in angular size and captures more light from the point.  
 The size of the image of the point is NOT changed.  
*Entrance pupil determines brightness*



Inverse square rule for point sources  
 Image illuminance  $\propto 1 / d^2$

**For extended sources (large objects)**  
 Object does not get brighter  
*Retinal image of object gets larger,  
 but average retinal illuminance remains the same*


### Bottom line on brightness



**FOR POINT SOURCES**  
 Image brightness determined by the area of light capture  
*Either TS objective OR M times eye pupil,  
 (whichever is smaller)*

**FOR EXTENDED SOURCES**  
 Image brightness determined by size of beam to enter eye  
*Either eye pupil OR exit pupil of TS,  
 (whichever is smaller)*

### All images from the system are within the Image Field Cone



FoV  
20 cm at 60 cm

Field limiter  
10 cm diameter  
Aperture in a card

FoV  
5 cm at 15 cm

Eye-to-aperture  
= 30 cm

IFAR = 10 cm/30 cm  
 = 0.33 = 19°

### FIELD ASPECT RATIOS


**Image Field Aspect Ratio (IFAR)**  
*Is the ratio of the diameter of image space aperture that limits the field to its distance from the eye point.*

**Object Field Aspect Ratio (OFAR)**  
*Is the ratio of the diameter of object space aperture that limits the field to its distance from the object space representation of the eye point*

*For Telescopes, IFAR = M (OFAR)*

### Field of View for simple Galilean and Keplerian telescopes

**Keplerian**  
 IFAR =  $A_2 / z$   
 or IFAR =  $A_2 / (d/M)$



**Keplerian**  
 Field limitation by Ocular lens (or field lens?)  
 Ocular lens subtends a large angle at the eye

**Galilean**  
 Field limitation by image of objective (exit pupil of TS)  
 TS exit pupil is small and remote. Thus subtends a small angle at the eye

**Galilean**  
 IFAR =  $(A_1 / M) / (z + d/M)$

### Rule of thumb for Fields of View of Telescopes


**Keplerian telescopes** (Usually IFAR = 1)  
 $FoV \approx EVD$

**Galilean Telescopes** (Usually IFAR =  $1/M_o$ )  
 $FoV \approx EVD / M_{ts}$

**Depth of Field of Telescopes**  
 $DoF \approx EVD^2$

4x12 TS with +5.00D cap  $EVD = (20 \text{ cm} / 4) = 5 \text{ cm} = 0.05 \text{ m}$

$DoF = (0.05)^2 = 0.0025 \text{ m} = 2.5 \text{ mm}$



### Imaging through afocal telescopes

*finite object distances*

Height of Image is **SMALLER**  
than Height of Object  
 $h' = h/M_{ts}$

Image Angle at Exit Pupil is **LARGER**  
than Object Angle at Entrance Pupil  
 $\omega' = M_{ts}\omega$

Image Smaller by  $M$   
and Closer by  $M^2$   
*(exact, if referenced to exit pupil)*

### Imaging through telescopes

Image Smaller by  $M_{ts}$ ..... and Closer by  $M_{ts}^2$   
*(referenced to exit pupil)*

**Image distance**  
relative to exit pupil  $v = u/M_{ts}^2$   
relative to ocular  $v = u/M_{ts}^2 + d/M_{ts}$

**Image Vergence**  
 $V_{xp} = M_{ts}^2 U_{obj}$  *(referenced to exit pupil)*

**Fried's formula** Emerging Vergence from Ocular  
 $V_{oc} = M_{ts}^2 U_{obj} / (1 - dM_{ts} U_{obj})$

$V_{oc} \approx M_{ts}^2 U_{obj}$

### Vergence amplification

**Implications**

Small change in object distance  
=> large change in accom demand

Focus range of TS might be insufficient to view close  
*May need to place chart at 4 m and use +0.25 D lens*

Weak lens on objective => large change in BVP

### Vergence amplification

**Implication**  
Risky to measure BVP of TS on lensometer


Lensometer for BVP of TS  
difficult to ensure lensometer is in precise adjustment  
*(zero vergence should enter the "lensometer eyepiece")*

Any error of lensometer-eyepiece focus is amplified by  $M_{ts}^2$  of TS

### Methods for checking telescopes


- 1 Direct Comparison (G or K)
- 2 Exit pupil measurement (G or K)
- 3 Vergence change at lensometer (G only)
- 4 Longitudinal magnification (K only)
- 5 Lateral Magnification -direct (G or K)
- 6 Lateral Magnification - lensometer (G or K)

**1 Comparison Method**




Telescope view

**2 Entrance/Exit Pupil Ratio**




**3. Lensometer method for measuring telescope magnification**

Galileian Only



**3. Lensometer method for measuring telescope magnification**

**3 Vergence change method ( Galileian )**




METHOD  
 Measure diff between BVP with and w/o + 1 D on obj lens  
 Estimate length of Galileian ts  
 Use table to determine magnification

**3. Lensometer method for measuring telescope magnification**

BVP when +1.00 D added to TS on Objective	MAGNIFICATION FOR Telescope Length in mm =					
	10	20	30	40	50	60
4.00	1.98	1.96	1.94	1.92	1.90	1.88
5.00	2.21	2.19	2.16	2.14	2.11	2.09
6.00	2.42	2.39	2.36	2.33	2.30	2.28
7.00	2.61	2.58	2.54	2.51	2.48	2.44
8.00	2.79	2.75	2.71	2.67	2.64	2.60
9.00	2.96	2.91	2.87	2.83	2.78	2.74
10.00	3.11	3.06	3.02	2.97	2.92	2.88
11.00	3.26	3.21	3.16	3.10	3.05	3.00
12.00	3.40	3.35	3.29	3.23	3.18	3.12
13.00	3.54	3.48	3.42	3.35	3.30	3.24
14.00	3.67	3.60	3.54	3.47	3.41	3.35
15.00	3.80	3.73	3.65	3.58	3.52	3.45
16.00	3.92	3.84	3.77	3.69	3.62	3.55
17.00	4.04	3.96	3.88	3.80	3.72	3.64
18.00	4.15	4.07	3.98	3.90	3.82	3.74
19.00	4.26	4.17	4.08	4.00	3.91	3.83
20.00	4.37	4.28	4.18	4.09	4.00	3.91


**4 Longitudinal Magnification (Keplerian)**



$d' = d M^2$   
 So  $M = \text{SQRT} (d'/d)$





Thin test plate close to ocular  
 Closer than exit pupil  
 Distinctive marks on both surfaces

**5 Lateral Magnification (Keplerian)**



$h' = M h$

Enlarged image of scale Located close to objective lens  
 Transparent scale located just inside exit pupil



- ## Dangle-ometry
- ### Determining Telescope Mounting Angle
1. Select telescope viewing point TVP (center of hole)  
*Must be at least half diameter plus 3mm from eyewire*
  2. Locate a fixation point at eye level  
*Ensures horizontal line of sight*
  3. Subject tilts head to align eye, TVP and fixation point
  4. Protractor with plumb line *measures spectacle plane angle*  
  
*When line of sight through TVP is horizontal,  
THEN Spectacle Frame Angle = Mounting Angle*

### OPTICAL PARAMETERS of a few telescopes

Manufacturer's Specifications	Brand	Magnification	Measured		Optical Factory			Optical Parameters		Physical		Size		Field	
			mm	mm	mm	mm	mm	mm	mm	mm	mm	mm	mm	mm	mm
	Walters	2.75x	2.9	0.60	2.7	14	>-20+11	812	3600	32	39	48	22	33	12
	Walters	4x12	3.9	0.87	2.9	22	>-20+20	1230	25300	47	47	58	30	47	12
	Walters	6x16	5.9	1.07	2.7	27	>-20+20	1612	17300	67	69	81	31	57	10
	Walters	7x25	6.8	1.20	3.0	49	>-20+5.5	2177	15100	157	113	136	42	62	9
	Walters	8x32	7.8	1.02	3.9	66	>-20+8	3021	13300	194	121	159	41	54	7
	Specewell	2.75x	2.7	0.67	3.1	13	>-20+12	813	37100	34	49	54	15	27	14
	Specewell	3x9	3.0	0.60	3.0	21	>-20+20	931	33301	32	44	50	24	33	11
	Specewell	4x12	3.9	0.87	2.8	19	>-20+20	1131	26300	38	46	58	29	47	12
	Specewell	6x16	5.7	1.10	2.7	22	>-20+20	1612	17500	67	89	129	32	54	7
	Specewell	10x20	9.5	1.00	2.0	38	>-20+20	1912	11100	86	105	136	32	53	6
	Specewell	15x20	9.8	1.10	3.0	43	>-20+20	3015	10100	172	148	176	47	57	6
	Sebi	6x18	6.0	0.91	2.2	800	>-20+20	1313	19117	125	67	68	32	50	49
	Sebi	6x15	5.6	0.83	2.6	201	-2.5+13	1512	20113	51	54	84	35	35	45
	Sebi	8x20	7.9	1.03	2.5	29	>-20+20	2012	13100	79	89	129	32	54	7
	Sebi	10x20	9.5	0.68	1.9	440	-15+20	197	11105	125	82	83	44	33	38
	Mayflower	6x15	6.0	0.83	2.5	234	>-20+20	1513	17305	54	57	59	37	45	8
	Mayflower	8x30	7.4	1.08	4.0	246	>-20+20	3015	14101	201	110	116	49	50	57
	Deigoen IV	2.0x EFES	2.0	0.72	6.8	18	-16+2.3	1416	50100	30	45	52	23	39	19
	DFV	3.0x EFES	2.8	0.79	4.6	23	-14+4	1316	36100	30	46	55	23	43	15
	DFV	4.0x EFES	4.2	0.77	4.8	27	>-20+4	2016	24100	42	59	81	26	42	10
	DFV	5.0x EFES	5.2	0.85	3.8	38	>-20+5	3016	19100	72	55	81	26	46	9
	DFV	6.0x EFES	6.9	0.84	2.9	31	>-20+2.5	2011	15100	35	55	82	25	45	7
	DFV	7.0x EFES	8.0	0.92	2.5	35	>-20+8	2010	13100	38	58	82	26	50	6
	DFV	8.0x EFES	8.1	0.91	2.5	35	>-20+7	2011	12601	35	67	82	25	49	6
	DFV	10x EFES	9.9	0.64	2.0	32	>-20+6.5	2011	10100	41	61	88	26	36	4
	Beecher	4x20	4.2	0.97	4.2	77	-13+11	1012	24002	116	54	57	41	47	55
	Beecher	5x25	5.6	1.05	5.1	175	-11+12	2912	18103	94	40	42	33	59	55
	Beecher	7x30	6.8	1.04	4.5	224	-9+10	3012	15101	105	38	43	35	59	55
	Beecher	10x35	9.2	0.99	3.8	423	-20+1.8	3512	11102	137	76	79	38	59	53